

Product highlights

- Effective infusion care by hanging safely IV-bottles and bags
- IV pole can be installed on either side of chair
- Very easy and fast assembly by a plug in and out system using one hand-screw only
- Space-saving infusion care (no need for any additional mobile IV pole stands)
- Comfortable height adjustment of bottle holder through a hand-screw height adjustment with an integrated, surface-gently mold clamping element
- All plastic parts are light-resisting
- Integrated pull-out guard avoids accidental extraction of the adjustable tube
- A sealing cap prevents the IV-pole from liquids getting inside

Technical data

Load capacity per hook	2 kg
Load capacity per holder	8 kg
Material	Stainless steel tube, plastic bottle holder
Type of height adjustment	stepless height adjustment by hand-screw
Arrangement of the hooks	90° to each other
Number of hooks	4 pcs
IV pole height adjustment	62 cm
Available for	Comfort-2 Eco/Flex Comfort-3 Eco/Flex Comfort-4 Eco/Battery/Plus Therapy chairs

Infusion bottle holder
with stepless height adjustment
installed on Comfort-2 Therapy Chair



Infusion bottle holder
on *Comfort-3 Therapy Chair*



Infusion bottle holder
on *Comfort-4 Therapy Chair*



